
UreSil® POSI-STICK® 2 Introducer Set with Radioopaque Marker

INTENDED USE:

This device is intended for interventional procedures: specifically a Seldinger technique placement of a percutaneous catheter.

CONTRAINDICATIONS:

This set is not intended for coronary or cerebrovascular use.

REF	PS2-38	PS2-38SW	PS2-35	PS2-35SW	PS2-00	PS2-18*
0.018" (4.6mm) GUIDEWIRE WITH FLOPPY RADIOPAQUE TIP	✓	✓	✓	✓	✓	✓
0.035" (0.89mm) / 0.038" (0.96mm) GUIDEWIRE WITH FLOPPY TIP & 9" TIP	✓	✓	✓	✓	✓	✓
0.035" (0.89mm) / 0.038" (0.96mm) SUPER-STIFF STRAIGHT GUIDEWIRE	✓	✓	✓	✓	✓	✓

*Replacement: 0.018" (4.6mm) Guidewire for PS2-35 & PS2-38

INSTRUCTIONS FOR USE:

1. Insert the 21 GA trocar-style needle using standard technique and guide into the drainage site.
2. Remove the trocar leaving the needle in place.
3. Advance the .018" (4.6 mm) guidewire through the needle into the drainage site.
4. Remove the needle while maintaining the .018" (.46 mm) guidewire in position.
5. Advance the "POSI-STICK" over the .018" (.46 mm) guidewire into the desired position.
6. Release and remove the stiffening cannula and the 4 Fr. dilator, leaving the 6 Fr. sheath and .018" (.46 mm) guidewire in position.
7. Advance the .038" or .035" (.96 mm or .89 mm) guidewire through the 6 Fr. sheath (along side the .018"/.46mm guidewire) into the drainage site.
8. Remove the .018" (.46 mm) guidewire.
9. Remove the 6 Fr. sheath.

CAUTION:

1. Withdrawal of guidewires should be smooth and without resistance. If resistance is felt, carefully remove wire(s) and system as a unit.
2. Avoid bending, kinking, and uses other than percutaneous catheter placement.
3. Guidewires should be routinely inspected prior to use. Should any deformities be present, discard. Do not use deformed guidewires.
4. Should resistance occur during insertion, DO NOT advance the guidewire.
5. Avoid withdrawing a guidewire through a metal needle. The teflon coating can shear against the needle bevel.

STERILE EO 2 CE 0086

Sterile if package is unopened and undamaged. Non-Pyrogenic. For single use only. The reuse of this single-use device can lead to patient infection and/or device malfunction.

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician. UreSil and POSI-STICK are registered trademarks of UreSil, LLC.

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