Centesis Catheter

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	5
Anwendungsgebiete	6
Ενδείξεις χρήσης	7
Felhasználási javallatok	8
Indicazioni per l'uso	9
Bruksindikasjoner	10
Indicações de utilização	11
Indicaciones de uso	12
Indikationer för användning	13
Kullanım Endikasyonları	14

Centesis Catheter

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION: The Centesis Catheter is an over-the needle catheter for percutaneous aspiration and small volume drainages for durations not to exceed 24 hours. Each Centesis Catheter includes the following components: Catheter with three distal sideholes, percutaneous entry needle and porous plug.

INTENDED USE

The Centesis Catheter is intended for use in percutaneous fluid aspirations and small volume drainage procedures.

CONTRAINDICATIONS

- Where percutaneous drainage catheterization is unacceptable
- Bleeding diathesis and uncontrolled hypertension
- Anticoagulant use
- Chemical ablation of tissue

WARNINGS

- The Centesis Catheter is supplied sterile for single use only.
- Do not exceed insertion dwell time of 24 hours.
- Do not alter this device in any way.
- Do not resterilize the device. Re-sterilization could change the physical characteristics of the materials and should not be attempted.
- Do not reuse the device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- During insertion avoid contact with bone, cartilage and scar tissue, which can damage the needle and catheter tip.

PRECAUTIONS

The Centesis Catheter deployment procedure should be performed by physicians thoroughly trained in drainage techniques and procedures.

Do not use the Centesis Catheter if the packaging is damaged.

Inspect the Centesis Catheter prior to use for any bends or kinks. Do not use a damaged needle because trauma to the patient and/or inability to place the drainage catheter may occur.

COMPLICATIONS: These are generally associated with drainage procedures and may occur when using the Centesis Catheter:

- Local or systemic infection
- Pneumothorax
- Inadvertent arterial or venous puncture resulting in bleeding complications
- Catheter or needle fracture and embolization

CLINICAL PROCEDURE

The Centesis Catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

PREPARATIONS FOR USE

1. Carefully inspect the Centesis Catheter packaging and components for damage prior to use.
2. Using sterile technique, open the pouch and transfer the Centesis Catheter into the sterile field.
3. Remove the needle from the catheter. Thoroughly flush both the needle and the catheter with saline, then reassemble them.
4. Prep skin and drape area of anticipated puncture as required.

DEPLOYMENT PROCEDURE

Deploy the Centesis Catheter according to the following steps:

1. Prepare the skin incision site using local anesthetic.
2. Introduce the Centesis Catheter into the cavity using ultrasound guidance, CT or fluoroscopy.
3. Remove the needle while leaving the catheter in position.
4. Connect a syringe or suction tube to the hub of the catheter to collect the drainage.
5. Upon completion of the drainage, remove the catheter and discard according to institutional practice.

WARNING: During insertion avoid contact with bone, cartilage and scar tissue, which can damage the needle and catheter tip.

PACKAGING & STORAGE

The Centesis Catheter has been sterilized with ethylene oxide.
Handle with care.

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Galt Medical Corp
U.S.A.
UreSil, LLC