
INTENDED USE:
 The THORA-VENT Thoracic Vent Procedure Tray has been designed for the treatment of spontaneous, traumatic or iatrogenic. It can also be used as a pleural catheter to utilize for evacuation of air from the pleural space after thoracic, or after treatment of pneumothorax which may result after chest tube removal. Do not use the drain for continuous use until approved by your physician. Intended for short-term use (more than 24 hours and less than 30 days).

TROCAR

INSTRUCTIONS FOR USE:

- Select a site, preferably at the second interspace in the mid-axillary line or at another appropriate location.
- Prepare the skin by scrubbing with alcohol. Using the scalpel, make a small skin incision at the selected site. If desired, the adhesive patch may be applied directly adjacent to the incision.
- Remove the trocar from the trocar and introduce the trocar into the device through the self-sealing port. The point of the trocar should extend just beyond the tip of the catheter.
- Remove the protective sheath from the catheter. Peel away the center portion of the paper cover from the adhesive patch. Do not retract the trocar unless the catheter has been completely introduced. The trocar/catheter assembly into the pleural space through the prepared incision staying immediately above the top border of the rib. Do not retract the rib or its entrapment.
- Entry of the catheter into the pleural cavity will be demonstrated by fluctuation of the red signal in the diaphragm. Then the tip of the catheter enters the pleural space, stop advancing the trocar and introduce the full length of the catheter. Remove the trocar. Do not retract the trocar unless the catheter has been completely removed from the patient. While the tip of the catheter is removed, the red signal diaphragm will deflect upwards when a pressure of greater than ($>$ 2 mm Hg (2.72 cm water)) is transmitted through the catheter and will retract when a pressure of less than ($<$ 2 mm Hg (2.72 cm water)) is transmitted through the catheter. The indicator will not properly function if the patient is ventilated or is on constant suction.
- Apply the paper on the side flaps of the adhesive patch and adhere the patch to the chest wall. The skin surface must be completely dry to ensure proper adhesion of the patch.
- When the self-sealing port is not in use (not connected) with the tethered cap.
- The catheter is designed to be active, sutures should be placed through the adhesive patch and the side holes on the top of the device to further anchor it to the patient's skin.
- Do not use the catheter for suction; the device, do not tape over the air vents (see diagram).

CHECK FOR CONTINUOUS AIR LEAK:
 The thoracic vent occlusion plug can be used to check for a continuous air leak from the patient's lung. To occlude the Thoracic Vent, insert and secure the occlusion plug to the luer-lock on the side of the Thoracic Vent. The occlusion plug will seal the system will not more than 0.02 cm of air at a pressure of 15 mm Hg. Warning: Remove the occlusion plug when the thoracic vent is ventilated. The occlusion plug will not function if the patient is on suction. Seal the self-sealing port with the tethered cap after the occlusion plug is removed.

ASPIRATION:
 A LF can be manually pumped out of the pleural space using the aspiration cannula. Insert the aspiration cannula into the vent through the self-sealing port. Attach the syringe to the aspiration cannula using the luer fitting, and utilize the syringe to evacuate air from the pleural space. Do not attach the syringe to the self-sealing port. Remove the self-sealing port using a heavy twisting and pulling on the flange (see diagram) of the aspiration cannula. Immediately seal the self-sealing port with a tethered cap after removing the aspiration cannula. If small amounts of naturally occurring fluid accumulate in the thoracic cavity, the fluid can be removed via the drainage port using a syringe.


EXTERNAL SUCTION ATTACHMENT:
 If external suction is desired, remove the suction tubing set from its pouch. Close the clamp on the tubing set while the set is being connected. Insert the cannula through the self-sealing port of the Thoracic Vent and secure the tubing by twisting the luer lock. Attach the funnel to the suction system and then open the clamp. Appropriate evacuating pressures should be used. Do not connect or disconnect the thoracic vent vacuum. Always clamp the suction tubing set when suction is not being applied to the Thoracic Vent. If the suction tubing set becomes disconnected, remove it and immediately seal the self-sealing port with the tethered cap.

WARNINGS/COMPLICATIONS:
 Do not use the device for aspiration. Introduce through or removal from the self-sealing port on the Thoracic Vent.
 Do not attempt to reinsert the trocar unless the catheter has been completely removed from the patient.
 DIFFICULT CATHETER REMOVAL: If the catheter is inadvertently introduced through the rib perforation on the trocar to the pleural space, the thoracic vent may be removed using the following procedure: If the catheter is placed on the chest wall, the catheter may be damaged and locked within the patient.
 Removal of the catheter should be resisted. If resistance is encountered, the catheter may be locked. If a patient complains of chest pain, it is a sign of a catheter dislodgement. In this case, it is recommended that a restrained catheter be released by direct surgical incision by cutting the catheter at the chest wall. Immediately remove the device, securing it and percutaneously passing clamps over the catheter until it is released from its attachment.

NOT FOR CONTINUOUS DRAINAGE:
 The THORA-VENT Thoracic Vent is not indicated for hemithorax or fluid evacuation other than air. Even though the device can accommodate small amounts of fluid (5 cc) without affecting function, large amounts of fluid can affect device function. When using insertion a hemithorax, or other liquid collection is present, alternate drainage is indicated.

MAINTAIN PROPER SEAL:
 An air-tight seal must be maintained to prevent subcutaneous emphysema.
 Continuous drainage of air over extended periods should alert a consideration for additional interventional treatment.
 The physician should be aware of complications associated with the treatment of pneumothorax including re-expansion and laceration of intercostal vessels.
 Always clamp the suction tubing set when suction is not being applied to the Thoracic Vent.
 Do not disinfect the Thoracic Vent with alcohol (propanol). Alcohol (propanol) degrades the device.
 If the thoracic vent plug was used to perform the air leak check, remove the occlusion plug when the check is complete. The thoracic vent will vent air when the occlusion plug is removed.

Dispose of sharps using approved sharps container in accordance with applicable regulatory and institutional policy.

Sterile if package is unopened and undamaged. The reuse of this device can lead to patient infection and/or device malfunction. UreSill and THORA-VENT are registered trademarks of UreSill, LLC. U.S. Patent no. 6,770,063.
 TROCACAR: CE 140935557 (04/13/05) only
 Fax: 1-847-962-8105
 www.uresill.com

Manufactured by: Skottell, LLC, 610077 U.S.A.
 5418 W. Touhy Ave., Skokie, IL, 60077 U.S.A.

Authorized Representative in European Community: EMERGO EUROPE Private Limited, 20, 2514 AF The Hague The Netherlands

Contains or presence of (phthalates) benzyl butyl methyl cellosolve (PBZ)

MR Conditional
 PHT
 BBP

Contains or presence of (phthalates) benzyl butyl methyl cellosolve (PBZ)

MR Conditional
 PHT
 BBP

Content of non-sterile (phthalates) benzyl and butyl cellosolve (PBZ)

MR Conditional
 PHT
 BBP

Content of non-sterile (phthalates) benzyl and butyl cellosolve (PBZ)

MR Conditional
 PHT
 BBP

Content of non-sterile (phthalates) benzyl and butyl cellosolve (PBZ)

MR Conditional
 PHT
 BBP

Content of non-sterile (phthalates) benzyl and butyl cellosolve (PBZ)

MR Conditional
 PHT
 BBP

Content of non-sterile (phthalates) benzyl and butyl cellosolve (PBZ)

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